FH ORTHOPEDICS

TIGE FEMORALE SANS CIMENT ESOP 2

ESOP 2 CEMENTLESS FEMORAL STEM

ZEMENTFREIER FEMURSCHAFT ESOP 2

VÁSTAGO FEMORAL SIN CEMENTO ESOP 2

STELO FEMORALE SENZA CEMENTO ESOP 2

HASTE FEMORAL SEM CIMENTO ESOP 2

ESOP 2 FEMORALE STEEL ZONDER CEMENT

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Instructions for use
ESOP 2 CEMENTLESS FEMORAL STEM

WARNINGS AND INDICATIONS FOR USE
PLEASE READ THIS DOCUMENT CAREFULLY

Please request help from the distributor's sales department if you need additional information such as the surgical techniques which show details about implantation. In any case, refer to the documents provided.

DESTINATION

ESOP 2 femoral stems are intended to be implanted in the femoral region during total hip arthroplasty.

DESCRIPTION

The ESOP 2 femoral stem is an anatomical stem with a cervicodiaphyseal angle of 130° that exists in both a left-sided (pink coloured) and right-sided (green coloured) version.

This cementless stem is coated with a T40 + HAP coating. It is intended to be used without cement.

The ESOP 2 stem is available in the following range:

Size	1	2	3	4	5	6	7
ESOP 2 right (green)	✓	✓	✓	✓	✓	✓	✓
ESOP 2 left (pink)	✓	✓	✓	✓	✓	✓	✓

FH INDUSTRIE FEMORAL HEADS

Taper 12/14	Necklength (mm)	-5	-3.5	+0	+2	+3.5	+7
	Diameter (mm)						
STAINLESS STEEL and CoCr	Ø22.2	×	×	✓	✓	×	×
	Ø28	✓ ⁽¹⁾	✓	✓	×	✓	✓
	Ø32	×	✓	✓	×	✓	✓
⁽¹⁾ stainless steel heads only							
BIOLOX® FORTE	Necklength (mm)	-5	-3.5	+0	+2	+3.5	+7
	Diameter (mm)						
	Ø28	×	✓	✓	×	✓	×
	Ø32	×	✓	✓	×	✓	×
	Ø36	×	✓	✓	×	✓	×

FEMORAL HEADS MANUFACTURED BY CERAMTEC

Taper 12/14	Necklength (mm)	-3.5	+0	+3.5	+7
BIOLOX® DELTA	Diameter (mm)				
	Ø28	28 12/14 S Ref: 38.49.7175.445.00	28 12/14 M Ref: 38.49.7175.455.00	28 12/14 L Ref: 38.49.7175.465.00	*
	Ø32	32 12/14 S Ref: 38.49.7175.665.00	32 12/14 M Ref: 38.49.7175.675.00	32 12/14 L Ref: 38.49.7175.685.00	32 12/14 XL Ref: 38.49.7181.345.00
Ø36	36 12/14 S Ref: 38.49.7179.275.00	36 12/14 M Ref: 38.49.7179.285.00	36 12/14 L Ref: 38.49.7179.295.00	*	

MATERIALS THE IMPLANTS ARE MADE FROM

The materials that the implants are made from are noted on each product label.

The implant is radio-visible.

COMPATIBILITY BETWEEN DEVICES

It is imperative to assemble:

- heads and stems having the same degree of taper (e.g., 12/14)
 - cemented (AC) liners or cups with femoral heads of the same diameter (for example: Ø28)
 - biologically fixed implants ("SC" meaning cementless), for which the type of biological fixation is specified on the product label, must not be cemented.
 - the femoral stems are compatible with BIOLOX® DELTA heads, for which the references are cited in the "CERAMTEC FEMORAL HEADS" paragraph.
 - the highly cross-linked TRIANON liner should only be fixed with an alumina head, preferably BIOLOX® FORTE or BIOLOX® DELTA; or a Cobalt chrome (CoCr) metal head manufactured exclusively by FH Industrie.
- Under no circumstances should a stainless steel head be implanted against a highly cross-linked TRIANON liner**

Materials compatibility

for heads/stems/cups

(excluding degree of taper):

	Cups			Stems
	AC cups or SC cup liners or AC cup liners			PHARO reinforced cup
	UHMWPE	TRIANON PE (highly cross-linked UHMWPE)	BIOLOX® DELTA Alumina	UHMWPE 316LVM
ESOP 2 stem				
STAINLESS STEEL heads	✓	×	×	✓
CoCr heads	✓	✓	×	✓
BIOLOX® DELTA or BIOLOX® FORTE alumina heads	✓	✓	✓	✓
Cephalic heads	×	×	×	✓

The ancillary device to be used for the insertion of hip implants is specified in the surgical technique.

PERFORMANCE

The device helps restore joint function and withstands the usual stresses to which the joint in question is submitted.

INDICATIONS FOR USE:

ESOP 2 femoral stems are indicated for the following aetiologies:

- Hip disorders (various forms of OA and RA of the hip, Trauma to the head or the acetabulum and its sequelae, AVN of the femoral head with or without osteoarthritis)
- Femoral neck fracture

CONTRAINDICATIONS

- 1 - Infection or latent infection
- 2 - A mental or neuromuscular disorder that might create an unacceptable risk of instability, prosthetic fixation failure, or complications after surgery
- 3 - Insufficient bone stock.
- 4 - Known allergy to any component of the material mentioned on the product label
- 5 - Metabolic diseases that might compromise bone regrowth
- 6 - Non-cooperative patient who is unable to follow recommendations

PRECAUTIONS

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.

The outcome of arthroplasty depends on the patient's medical history. It is indispensable that the patient be psychologically prepared.

Patients receiving a hip joint prosthesis must be informed of the limitations of the prosthesis, including among others, the effect of their weight and activity level on the lifespan of the prosthesis. They must be advised as to how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device location.

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated all aspects of the surgery, as well as the limitations of the device.

- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.

- The contact zones of the implant must be clean and free of debris (bone, cement) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.

- An alumina head cannot be revised with another alumina head unless the corresponding stem is revised as well.

- Care should be taken to ensure that the presence of other devices will not harm the functional integrity of the device.

- Do not combine implants from various origins. The only possible combinations of components are indicated in the section

COMPATIBILITY BETWEEN DEVICES.

- If important information on the label (marketing ref., size, lot number) and/or if the laser marking (lot number, taper diameter, manufacturer) on the implant is illegible, do not use the implant.

CEMENTLESS IMPLANTS

During handling, great care should be given to the coating. Any contact with the coating before implantation must be avoided. Biologically fixed implants must not be cemented. The type of coating is identified on the label.

IMPLANTATION TECHNIQUE

The specific ancillary equipment, delivered non-sterile, is necessary for implanting implants designed by FH Industrie and distributed by FH ORTHOPEDICS or its distributors. The surgical techniques are available from FR ORTHOPEDICS or its distributors. FH ORTHOPEDICS or its distributors can supply radiological templates that the surgeon can use to determine the implant size.

CONDITIONS OF USE

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined above.

WARNINGS

- Patients receiving a prosthesis must be informed that the lifespan of the implant may depend on their weight and their activity level.
- Damaged or defective implants must not be used.
- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics.
- Using an implant that is of the wrong size or lateralization runs the risk of decreasing its stress resistance.
- The components are part of a modular system. FH Industrie is not responsible for use of another manufacturer's components, (apart from femoral heads from the manufacturer CERAMTEC for which the references are noted in the table "FEMORAL HEADS MANUFACTURED BY CERAMTEC"), or use of components for purposes other than those intended.
- The implants must not be modified or subjected to any treatment.
- No chemical products (ex. alcohol) should be applied to the device during clinical practice.

ADVERSE EFFECTS

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following risks:

- pain
- loss of mobility
- wear of the UHMWPE components. Wear particles, especially from components made of UHMWPE, may cause osteolysis that may require another surgery.
- implant rupture as a result of inappropriate activities, trauma, or other stress specific to the patient's activities.
- loosening or dislocation of an implant following inadequate initial fixation, latent infection, early or excessive stress, poor positioning of the components, or trauma.
- allergy to one of the compounds of the material noted on the product label.
- possible evolving, sometimes asymptomatic, bone resorption that may appear around prosthetic components as a result of reactions to a foreign body.
- hematoma, slowed healing, deep vein thrombosis, pulmonary thrombosis, vascular lesions
- pulmonary embolism
- infection, and possibly death
- a cardiovascular problem

These adverse effects may lead to another operation or a revision.

IMPORTANT INFORMATION TO GIVE THE SURGEON AND THE PATIENT:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

If undergoing an MRI exam, the patient must inform the staff that he/she has a metal implant:

- Metal implants carry a risk of migration and/or heating during an MRI exam.
- Metal implants can cause artifacts (signal distortion and/or signal loss) which may interfere with the interpretation of the imaging results.

PACKAGING AND STERILIZATION

The information on the product label ensures manufacturing traceability.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly closed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. In this case, the product must be returned. The product must not be re-sterilised by any method whatsoever at the risk of altering its mechanical characteristics. The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation (R). This colour may be altered by poor storage conditions: heat, humidity, light etc. In all cases, an indicator showing the pre-sterilisation colour, yellow, may indicate that the product is not sterile and, in this case, the product must be returned. Desterilised implants cannot be returned. The expiry date is noted on the product label.

DEVICE DISPOSAL

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE OF FIRST CE MARKING

ESOP 2 CEMENTLESS femoral stem: 2016

The reference text is the French text.

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UK, FH ORTHO LTD