GROUP

FH
ORTHO™

FH ORTHOPEDICS

PROTHESE DE CHEVILLE EASYMOVE

EASYMOVE ANKLE PROSTHESIS

EASYMOVE SPRUNGGELENKPROTHESE

PRÓTESIS DE TOBILLO EASYMOVE

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INSTRUCTIONS FOR USE

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Consult
instructions
for use

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**WARNINGS AND INDICATIONS FOR USE
PLEASE READ THIS DOCUMENT CAREFULLY**

Do not hesitate to consult your supplier's sales department if you need further information. This document is not exhaustive, nor is it a technical user manual containing detailed layout drawings. In any case, refer to the documents provided.

DESCRIPTION

FH ORTHOPEDICS provides an ankle prosthesis.

It is a semi-constrained prosthesis consisting of two complementary components: a metal tibial base associated with a polyethylene insert and a metallic resurfacing talar component.

X-ray templates are provided for selecting the size of the implant before surgery.

Instruments are also provided:

- Test parts for testing the implantation during surgery
- Instruments for the assembly and correct implantation of the prosthesis.

For a more detailed description of the implants and their use, refer to the technical documentation or contact your FH Orthopedics representative. It is essential to implant the total ankle prosthesis with the FH Orthopedics instrumentation that has been specifically designed for this purpose. FH Orthopedics implants must be assembled using FH Orthopedics components that are defined as being compatible with one another. Appropriate implants can be selected using the recommendations in the surgical technique, the test parts, and the templates provided with the instrumentation.

Symbols may be used to identify some implants (labelling or marking) with the following meanings: **L** = left; **R** = right; **TH** = thickness.


INDICATIONS FOR USE:

The total ankle prosthesis is indicated as a total ankle replacement for primary surgery in patients with ankle joints damaged by severe rheumatoid, post-traumatic or degenerative arthritis.

The components are intended for uncemented use.

COMPATIBILITY:

Recommendations on how to combine the components are noted in the following drawing:

PERFORMANCE:

The total ankle prosthesis should decrease pain and restore ankle function compared to the pre-operative state. This medical device is used for primary treatment.

MATERIALS THE IMPLANTS ARE MADE FROM

The materials that the implants are made from are noted on each product label.

The implants are radiopaque except for the UHMWPE insert (ultra-high molecular weight polyethylene).

The tibial base is made of titanium alloy (Ti6Al4V) in accordance with ISO 5832-3. The talar component is made of cobalt-chrome alloy (CoCr) in accordance with ISO 5832-4. The polyethylene insert is made of UHMWPE in accordance with ISO 5834-2. The tibial and talar implants are coated with porous titanium by plasma spray according to ASTM F1580.

CONTRAINDICATIONS:

- Infection
- Infection sequelae
- Systemic infection, fever and/or local inflammation
- Complete talar necrosis
- Insufficient bone stock or poor skin coverage around the ankle joint that would render the procedure unjustifiable
- Persistent cutaneous lesion

- Significant ligamentous laxity
- Severe osteoporosis
- Ankle arthrodesis with malleolar exeresis
- Neuromuscular or mental disorders which might jeopardise fixation and postoperative care
- Neurobiological disorders
- Nonfunctional lower limb muscles
- Complete loss of collateral ligament of the ankle
- Charcot arthropathy
- Distant foci of infection from genitourinary, pulmonary, skin and other sites, dental focus infection which may cause hematogenous spread to the implant site
- Bone immaturity
- Known allergy to one of the materials
- Vascular insufficiency at the ankle joint
- Inability of the patient to follow the surgeon's recommendations and the physical therapy program
- Patient pregnancy

PRECAUTIONS:

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.
- The outcome of ankle prosthesis depends on the patient's medical history. It is indispensable that the patient be psychologically prepared.
- Patients must be informed of the limitations of the device, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised as to how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint.
- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.
- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.
- The contact zones of the implant must be clean and free of debris (bone, cement) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.
- The interior of the tibial implant must be perfectly clean and free of all tissue or bony debris that might hinder the setting of the insert against it.

- It is advisable to make sure that the presence of other devices does not compromise the device's functional integrity.
- Implants provided by FH Orthopedics must not be mixed with implants of another origin, or with different systems.
- During handling, great care should be given to the coating. Any contact with the coating before implantation must be avoided. Biologically fixed implants must not be cemented. The type of coating is identified on the label.

IMPLANTATION TECHNIQUE:

The specific ancillary equipment, delivered non-sterile, is necessary for implanting implants designed by FH Industrie and distributed by FH ORTHOPEDICS or its distributors. The surgical techniques are available from FH ORTHOPEDICS or its distributors. For certain implants, FH ORTHOPEDICS or its distributors may provide x-ray templates to enable the surgeon to evaluate the size of the implant.

CONDITIONS OF USE

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by an orthopaedic surgeon who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined above.

WARNINGS:

- Patients receiving an implant must be informed that the implant's lifespan may depend on their weight and activity level.
- Damaged or defective implants must not be used.
- The design of this device does not allow for its reuse, at the risk of altering its mechanical characteristics.
- Using an implant that is of the wrong size or lateralisation runs the risk of decreasing its stress resistance.
- FH Industrie is not responsible for use of another manufacturer's components, or use of components for purposes other than those intended.
- The implants must not be modified or subjected to any treatment.

ADVERSE EFFECTS:

The surgeon is responsible for providing the patient with complete information before the surgical procedure, and in particular for informing them of the following risks:

- Complications related to any surgical intervention

- Rupture of one or several implants following inappropriate activities, trauma or other loads related to the patient's activities.
- Disassembly of one or several implants due to inadequate initial fixation, latent infection, early or excessive stress, poor positioning of the components, or trauma.
- Potential for progressive asymptomatic bone resorption around prosthetic components due to a foreign body reaction.
- Rupture of the prosthesis following an undetected production fault
- Wear of the polyethylene components (UHMWPE). Wear particles, especially from components made of UHMWPE, may cause osteolysis that may require another surgery.
- Side effects that could require a new surgical procedure or revision.
- Dislocation of the prosthesis
- Sensitivity or allergy to one of the compounds of the material listed on the product label.
- Infection.
- Slow or poor healing.
- Instability.
- Secondary talar necrosis.
- Neuropathy.
- Pain
- Haematoma, deep vein thrombosis, pulmonary embolism, cardiovascular problem
- Rarely death

IMPORTANT INFORMATION TO GIVE THE SURGEON AND THE PATIENT:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

If undergoing an MRI exam, the patient must inform the staff that he/she has a metal implant:

- Metal implants carry a risk of migration and/or heating during an MRI exam.
- Metal implants can cause artifacts (signal distortion and/or signal loss) which may interfere with the interpretation of the imaging results.

PACKAGING AND STERILISATION

The information on the product label ensures manufacturing traceability and indicates which sterilisation method has been used (if applicable).

The tibial and talar implants are supplied having been sterilised by gamma radiation and the insert by ethylene oxide. The expiry date for sterilisation and integrity of the packaging must be verified.

An implant whose packaging has been opened or damaged or whose expiration date has passed must not be used.

Every precaution must be taken to ensure sterility when opening the packaging of the implant and during implantation.

The instruments can be supplied sterile.

See the instructions about instruments for how to handle and sterilise non sterile instruments. The templates are supplied non sterile and should not be sterilised.

For all other information about instruments, see the instructions provided for that purpose.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches or blister packs and films) are perfectly sealed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. The product must not be resterilised by any method whatsoever at the risk of altering its mechanical characteristics.

The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation or green for EtO gas sterilisation (EtO). This colour may be altered by poor storage conditions such as heat, humidity, light, etc.

In all cases, an indicator showing one of the pre-sterilisation colours (yellow or purple) may indicate that the product is not sterile and, in this case, the product must be returned.

Desterilised instruments cannot be returned.

The expiry date is noted in the product label.

DEVICE DISPOSAL

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE OF FIRST CE MARKING

Date of first CE marking: 2019.

The reference text is the French text.

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UK, FH ORTHO LTD