# ACTIFIT™

Substitut méniscal biodégradable

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	Référence	Année
1	<b>Tissue ingrowth after implantation of a novel, biodegradable polyurethane scaffold for treatment of partial meniscal lesions.</b> Verdonk R, Verdonk P, Huyse W, Forsyth R, Heinrichs EL. Am J Sports Med. 2011 Apr;39(4):774-82.	2011
2	<b>The safety and short-term efficacy of a novel polyurethane meniscal scaffold for the treatment of segmental medial meniscus deficiency.</b> Efe T, Getgood A, Schofer MD, Fuchs-Winkelmann S, Mann D, Paletta JR, Heyse TJ. Knee Surg Sports Traumatol Arthrosc. 2012 Sep;20(9):1822-30.	2012
3	<b>Meniscal scaffolds: early experience and review of the literature.</b> Spencer SJ, Saithna A, Carmont MR, Dhillon MS, Thompson P, Spalding T. Knee. 2012 Dec;19(6):760-5.	2012
4	<b>Successful treatment of painful irreparable partial meniscal defects with a polyurethane scaffold: two-year safety and clinical outcomes.</b> Verdonk P, Beaufils P, Bellemans J, Djian P, Heinrichs EL, Huyse W, Laprell H, Siebold R, Verdonk R. Actifit Study Group. Am J Sports Med. 2012 Apr;40(4):844-53.	2012
5	<b>Two-year follow-up study on clinical and radiological outcomes of polyurethane meniscal scaffolds.</b> De Coninck T, Huyse W, Willemot L, Verdonk R, Verstraete K, Verdonk P. Am J Sports Med. 2013 Jan;41(1):64-72.	2013
6	<b>Polyurethane scaffold for the treatment of partial meniscal tears. Clinical results with a minimum two-year follow-up.</b> Bulgheroni P, Bulgheroni E, Regazzola C, Mazzola C. Joints. 2013 Mar 21;1(4):161-6. eCollection 2013 Oct-Dec.	2013
7	<b>Biodegradable polyurethane meniscal scaffold for isolated partial lesions or as combined procedure for knees with multiple comorbidities: clinical results at 2 years.</b> Kon E, Filardo G, Zaffagnini S, Di Martino A, Di Matteo B, Marcheggiani Muccioli GM, Busacca M, Maracci M. Knee Surg Sports Traumatol Arthrosc. 2014 Jan;22(1):128-34.	2014
8	<b>Polyurethane scaffold in lateral meniscus segmental defects: clinical outcomes at 24 months follow-up.</b> Bouyarmane H, Beaufils P, Pujol N, Bellemans J, Roberts S, Spalding T, Zaffagnini S, Maracci M, Verdonk P, Womack M, Verdonk R. Orthop Traumatol Surg Res. 2014 Feb;100(1):153-7.	2014
9	<b>Actifit synthetic meniscal substitute: experience with 18 patients in Brest, France.</b> Baynat C, Andro C, Vincent JP, Schiele P, Buisson P, Dubrana F, Gunepin FX. Orthop Traumatol Surg Res. 2014 Dec;100(8 Suppl):S385-9.	2014
10	<b>Cell distribution and regenerative activity following meniscus replacement.</b> Moran CJ, Atmaca S, Declercq HA, Cornelissen MJ, Verdonk PC. Int Orthop. 2014 Sep;38(9):1937-44.	2014
11	<b>Improvement in outcomes after implantation of a novel polyurethane meniscal scaffold for the treatment of medial meniscus deficiency.</b> Schuttler KF, Pylvättgen S, Getgood A, Rominger MB, Fuchs-Winkelmann S, Roessler PP, Ziring E, Efe T. Knee Surg Sports Traumatol Arthrosc. 2015 Jul;23(7):1929-35.	2015
12	<b>Partial meniscus substitution with a polyurethane scaffold does not improve outcome after an open-wedge high tibial osteotomy.</b> Gelber PE, Isart A, Erquicia JI, Pelfort X, Tey-Pons M, Monllau JC. Knee Surg Sports Traumatol Arthrosc. 2015 Jan;23(1):334-9.	2015
13	<b>Results of polyurethane implant for persistent knee pain after partial meniscectomy with a minimum of two years follow-up.</b> Martín-Hernández C, Ranera-García M, Díaz-Martínez JV, Muniesa-Herrero MP, Floría-Arnal LJ, Osca-Guadalajara M, García-Aguilera D. Rev Esp Cir Ortop Traumatol. 2015 Jan-Feb;59(1):44-51.	2015
14	<b>Actifit® scaffold implantation: Influence of preoperative meniscal extrusion on morphological and clinical outcomes.</b> Faivre B, Bouyarmane H, Lonjon G, Boisrenoult P, Pujol N, Beaufils P. Orthop Traumatol Surg Res. 2015 Oct;101(6):703-8.	2015
15	<b>The magnetic resonance aspect of a polyurethane meniscal scaffold is worse in advanced cartilage defects without deterioration of clinical outcomes after a minimum two-year follow-up.</b> Gelber PE, Petrica AM, Isart A, Mari-Molina R, Monllau JC. Knee. 2015 Oct;22(5):389-94.	2015
16	<b>Midterm follow-up after implantation of a polyurethane meniscal scaffold for segmental medial meniscus loss: maintenance of good clinical and MRI outcome.</b> Schüttler KF, Haberhauer F, Gesslein M, Heyse TJ, Figiel J, Lorbach O, Efe T, Roessler PP. Knee Surg Sports Traumatol Arthrosc. 2016 May;24(5):1478-84.	2016
17	<b>Association of a hi-tech with a bio-tech technique in the treatment of early osteoarthritis of the knee: a case report.</b> Uboldi FM, Ferrua P, Parente A, Pasqualotto S, Berruto M. Joints. 2016 Jan 31; 3(4): 221-224	2016
18	<b>Comparative Study of Collagen versus Synthetic-Based Meniscal Scaffolds in Treating Meniscal Deficiency in Young Active Population.</b> Bulgheroni E, Grassi A, Campagnolo M, Bulgheroni P, Mudhigere A, Gobbi A. Cartilage. 2016 Jan;7(1):29-38.	2016
19	<b>Treatment of Painful, Irreparable Partial Meniscal Defects With a Polyurethane Scaffold: Midterm Clinical Outcomes and Survival Analysis.</b> Dhollander A, Verdonk P, Verdonk R. Am J Sports Med. 2016 Oct;44(10):2615-2621.	2016
20	<b>Polyurethane-based cell-free scaffold for the treatment of painful partial meniscus loss.</b> Filardo G, Kon E, Perdissia F, Sessa A, Di Martino A, Busacca M, Zaffagnini S, Maracci M. Knee Surg Sports Traumatol Arthrosc. 2017 Feb;25(2):459-467.	2017
21	<b>Actifit polyurethane meniscal scaffold: MRI and functional outcomes after a minimum follow-up of 5-years.</b> Leroy A, Beaufils P, Faivre B, Steltzen C, Boisrenoult P, Pujol N. Orthop Traumatol Surg Res. 2017 Jun;103(4):609-614	2017
22	<b>Magnetic Resonance Imaging and Functional Outcomes After a Polyurethane Meniscal Scaffold Implantation: Minimum 5-Year Follow-up.</b> Monllau JC, Poggioli F, Erquicia J, Ramíñez E, Pelfort X, Gelber P, Torres-Claramunt R. Arthroscopy. 2018 May;34(5):1621-1627.	2018
23	<b>First Clinical Application of Polyurethane Meniscal Scaffolds with Mesenchymal Stem Cells and Assessment of Cartilage Quality with T2 Mapping at 12 Months.</b> Olivios-Meza A, Pérez Jiménez FJ, Granados-Montiel J, Landa-Solis C, Cortés González S, Jiménez Aroche CA, Valdez Chávez M, Renán León S, Gomez-García R, Martínez-López V, Ortega-Sánchez C, Parra-Cid C, Velasquillo Martínez C, Ibarra C. Cartilage. 2019 Aug 7:1947603519852415.	2019
24	<b>Polyurethane scaffold implants for partial meniscus lesions: delayed intervention leads to an inferior outcome.</b> Condello V, Dei Giudici L, Perdissia F, Screpis DU, Guerriero M, Filardo G, Zorzi C. Knee Surg Sports Traumatol Arthrosc. 2019 Oct 25.	2019
25	<b>Medial Meniscus Scaffold Implantation in Combination with Concentrated Bone Marrow Aspirate Injection: Minimum 3-Year Follow-up.</b> Akkaya M, Şimşek ME, Gürsoy S, Çay N, Bozkurt M. J Knee Surg. 2020 Feb 10.	2020
26	<b>Polyurethane Meniscal Scaffold for the Treatment of Partial Meniscal Deficiency: 5-Year Follow-up Outcomes: A European Multicentric Study.</b> Toanen C, Dhollander A, Bulgheroni P, Filardo G, Zaffagnini S, Spalding T, Monllau JC, Gelber P, Verdonk R, Beaufils P, Pujol N, Bulgheroni E, Asplin L, Verdonk P. Am J Sports Med. 2020;48(6):1347-1355.	2020

DISTRIBUTEURS  
DISTRIBUTORS

GRUPE  
**FH ORTHO**

**FR, FH ORTHO SAS**  
3 rue de la Forêt - Zone Industrielle  
BP 50009  
68990 Heimsbrunn CEDEX - FRANCE  
Tél. +33 (0)3 89 81 90 92  
Fax: +33 (0)3 89 81 80 11  
info@fhortho.com  
www.fhortho.com

**USA, FH ORTHOPEDICS INC.**  
OrthoEx  
7327 E Tierra Buena Lane  
Scottsdale, Arizona 85260 - USA  
Phone: +1 (412) 965-0950  
customerservice@fhortho-us.com  
www.fhortho.com

**PL, FH ORTHO POLSKA**  
Ul. Garbary 95/A6,  
61-757 Poznan - POLSKA  
Phone: +48 61 863 81 27  
Fax: +48 61 863 81 28  
biuro@implants24.pl  
www.fhortho.com