


FOOT

tictac

TOE INTERPHALANGEAL CAGE
TOE ARTHRODESIS CAGE



TICTAC $\varnothing 2,5 - 0^\circ$
ref. 270 390



TICTAC $\varnothing 2,8 - 0^\circ$
ref. 268 612



TICTAC $\varnothing 3,1 - 0^\circ$
ref. 268 614



TICTAC $\varnothing 3,4 - 0^\circ$
ref. 268 616



TICTAC $\varnothing 2,5 - 10^\circ$
ref. 270 391



TICTAC $\varnothing 2,8 - 10^\circ$
ref. 268 613



TICTAC $\varnothing 3,1 - 10^\circ$
ref. 268 615



TICTAC $\varnothing 3,4 - 10^\circ$
ref. 268 617



TICTAC $\varnothing 2,5 - 20^\circ$
ref. 270 392



TICTAC $\varnothing 2,8 - 20^\circ$
ref. 270 393



TICTAC $\varnothing 3,1 - 20^\circ$
ref. 270 394



TICTAC $\varnothing 3,4 - 20^\circ$
ref. 270 395

- > **Designed to fit all the patient's anatomy:**
 - 3 angulations: $0^\circ / 10^\circ / 20^\circ$
 - 4 diameters: 2,5mm, 2,8mm, 3,1mm and 3,4mm

- > **Intra-medullar press-fit** achieved by the implant design.

- > **Made of PEEK** as a way to improve the biocompatibility of load-bearing implants.

- > **Dedicated instruments**

> Instrumentation


DRILL Ø2

ref. 266 818




IPP CAGE TICTAC SHAPER P1 Ø2,5 BEPOD

ref. 271 066




IPP CAGE TICTAC SHAPER P2 Ø2,5 BEPOD

ref. 271 067




IPP CAGE TICTAC SHAPER P1 Ø2,8 BEPOD

ref. 268 618




IPP CAGE TICTAC SHAPER P2 Ø2,8 BEPOD

ref. 268 619



IPP CAGE TICTAC SHAPER P1 Ø3,1 BEPOD

ref. 268 620




IPP CAGE TICTAC SHAPER P2 Ø3,1 BEPOD

ref. 268 621



IPP CAGE TICTAC SHAPER P1 Ø3,4 BEPOD

ref. 268 622




IPP CAGE TICTAC SHAPER P2 Ø3,4 BEPOD

ref. 268 623




IPP CAGE TICTAC DRIVER HOLDER BEPOD

ref. 268 624



HOLDER BEPOD

ref. 271 147



DISTRIBUTEURS
DISTRIBUTORS

GRUPE
FH ORTHO™

FR, FH ORTHO SAS
3 rue de la Forêt - Zone Industrielle
BP 50009
68990 Heimsbrunn CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fhortho.com
www.fhortho.com

USA, FH ORTHOPEDICS INC.
OrthoEx
7327 E Tierra Buena Lane
Scottsdale, Arizona 85260 - USA
Phone: +1 (412) 965-0950
customerservice@fhortho-us.com
www.fhortho.com

PL, FH ORTHO POLSKA
Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fhortho.com



FABRICANT
MANUFACTURER

FR, FH INDUSTRIE
6 rue Nobel, Z.I. de Kernevez
29000 QUIMPER - FRANCE
Tél. +33 (0)2 98 55 68 95
Fax : +33 (0)2 98 53 42 13
contact-fhi@fhortho.com
www.fhortho.com