

TECHNICAL SHEET

Ext Liftfix

OUT IN - KNEE

ADJUSTABLE CORTICAL FIXATION SYSTEM
FOR ACL RECONSTRUCTION



Product delivered in double packaging
and sterilized by Oxide Ethylene.
Do not re-sterilize.

Product description

The device is composed of:

- A titanium plate (12 x 4 mm TA6V ELI),
- 1 adjustable loop (UHMWPE) - minimal length 10 mm),
- 2 pulling threads (1 blue/white and 1 white),
- 1 additional Jumbo titanium plate (TA6V ELI) for a larger contact surface with the cortical bone.

Surgical approach

- Compatible with standard surgical techniques used for ACL reconstruction with cortical fixation.



Plate:



Length: 12 mm
Width: 4 mm
Thickness: 1,3 mm
Tunnel width: 4,5 mm

Jumbo plate (optional):


Length: 20 mm
Width: 5,5 mm
Thickness: 2,6 mm
Tunnel width greater than
or equal to 5 mm









Pre-installation of the Jumbo plate (option 1)



Remove the pulling threads from the plate and reposition them on the Jumbo plate to allow its installation.



Placement of the Jumbo plate before tightening the implant (option 2)



Ref.	Designation
TM00140100	EXT LIFTFIX BUTTON TI
TM00130100	JUMBO BUTTON TI

GRUPE
FH ORTHO



DISTRIBUTEURS
DISTRIBUTORS

FR, FH ORTHO SAS
3 rue de la Forêt - Zone Industrielle
BP 50009
68990 Heimsbrunn CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fhortho.com
www.fhortho.com

USA, FH ORTHOPEDICS INC.
OrthoEx
7327 E Tierra Buena Lane
Scottsdale, Arizona 85260 - USA
Phone: +1 (412) 965-0950
customerservice@fhortho-us.com
www.fhortho.com

PL, FH ORTHO POLSKA
Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fhortho.com

Manufacturer:
Tulpar Medical solutions

CE 2292