



Technical sheet  
**REVISION OR ABLATION**


**ARROW I**

GROUPE  
**FH ORTHO**<sup>™</sup>

# CASE n°1

Replacement of ARROW I  
intermediate implants (identical)

Keep the stem  
**in place**  
**ARROW I**



Keep the baseplate or  
cemented glenoid **in place**  
**ARROW I**



The **ARROW I revision** ancillary **A266119** allows replacement of intermediate parts (insert and glenosphere)

**TO ORDER**



**ARROW II** ancillary + **ARROW I Revision Ablation** ancillary and associated implant nomenclatures



\* order universal anatomical ancillary ARROW I + the associated implant nomenclature

# CASE n°2

**Change** the stem and humeral intermediate implants (humeral insert or humeral head) **in ARROW II**



**ARROW II**



Keep the baseplate or cemented glenoid **in place**

**ARROW I**



**Replacement** of glenoid intermediate implants (insert or glenosphere) in

**ARROW I**



**TO ORDER**




**ARROW II** ancillary + **ARROW I Revision Ablation** ancillary and associated implant nomenclatures

\* order universal anatomical ancillary ARROW I + the associated implant nomenclature



# CASE n°3

**Change** of baseplate and glenoid intermediate implants (glenosphere or glenoid insert) **in ARROW II**


Keep the stem  
**in place**  
**ARROW I**



**Replacement** of humeral intermediate implant (humeral insert or humeral head) in **ARROW I**



**ARROW II**




**TO ORDER**



**ARROW II** ancillary + **ARROW I Revision Ablation** ancillary and associated implant nomenclatures

# CASE n°4

Change the whole platform in **ARROW II**



**ARROW I**

**ARROW II**

**TO ORDER**



**ARROW II** ancillary + **ARROW I Revision Ablation** ancillary  
and associated implant nomenclatures



# ARROW I IMPLANTS



## ARROW I REVISION ABLATION ANCILLARY (ref. A 266 119)

## ARROW I Revision implant nomenclature (ref. I 272 198)



### ARROW I HUMERAL INSERT

Ref.	Designation
257 057	ARROW humeral insert D36 H00
257 058	ARROW humeral insert D36 H05
257 059	ARROW humeral insert D36 H10
257 060	ARROW humeral insert D39 H00
257 061	ARROW humeral insert D39 H05
257 062	ARROW humeral insert D39 H10
257 063	ARROW humeral insert D42 H00
257 064	ARROW humeral insert D42 H05
257 065	ARROW humeral insert D42 H10



### ARROW I HUMERAL HEAD

Ref.	Designation
260 532	ARROW centered humeral head D48 H18
260 535	ARROW centered humeral head D50 H19
260 550	ARROW centered humeral head D54 H19




### ARROW I GLENSPHERE

Ref.	Designation
260 566	ARROW glenosphere D36
260 567	ARROW glenosphere D39
260 568	ARROW glenosphere D42

# ANATOMICAL METAL-BACK PROSTHESIS REVISION ARROW I

## ARROW I ANATOMICAL UNIVERSAL BASKET (ref. A 272 069)

## ARROW I implant nomenclature Anatomical revision (ref. I 272 104)



### CENTERED HUMERAL HEAD

Ref.	Diameter	Height (mm)
260 537	40	15
260 538	40	17
260 539	44	16
260 540	44	18
260 541	46	16
260 542	46	18
260 543	46	21
260 544	48	16
260 545	48	18
260 546	48	21
260 547	50	17
260 548	50	19
260 549	50	21
260 550	54	19
260 551	54	21



### OFF-CENTERED HUMERAL HEAD

Ref.	Diameter	Height (mm)
260 526	44	16
260 527	44	18
260 528	46	16
260 529	46	18
260 530	46	21
260 531	48	16
260 532	48	18
260 533	48	21
260 534	50	17
260 535	50	19
260 536	50	21



### GLENOID INSERT

Ref.	Size
260 556	44
260 557	46
260 558	48



DISTRIBUTEURS  
DISTRIBUTORS

GRUPE  
**FH ORTHO**<sup>TM</sup>

**FR, FH ORTHO SAS**  
3 rue de la Forêt - Zone Industrielle  
BP 50009  
68990 Heimsbrunn CEDEX - FRANCE  
Tél. +33 (0)3 89 81 90 92  
Fax : +33 (0)3 89 81 80 11  
info@fortho.com  
www.fortho.com

**USA, FH ORTHO INC.**  
OrthoEx  
7327 E Tierra Buena Lane  
Scottsdale, Arizona 85260 - USA  
Phone: +1 (412) 965-0950  
customerservice@fortho-us.com  
www.fortho.com

**UK, FH ORTHO LTD**  
Unit 1b, Century Park, Valley way  
Swansea Enterprise Park  
Swansea, SA6 8RP - UK  
Phone: +44 (0) 1792 464792  
Fax: +44 (0) 844 412 7674  
customer-servicesUK@fortho.com  
www.fortho.com

**PL, FH ORTHO POLSKA**  
Ul. Garbary 95/A6,  
61-757 Poznan - POLSKA  
Phone: +48 61 863 81 27  
Fax: +48 61 863 81 28  
biuro@implants24.pl  
www.fortho.com

**FR, FH INDUSTRIE**  
6 rue Nobel, Z.I. de Kernevez  
29000 QUIMPER - FRANCE  
Tél. +33 (0)2 98 55 68 95  
Fax : +33 (0)2 98 53 42 13  
contact-fhi@fortho.com  
www.fortho.com



FABRICANT  
MANUFACTURER