






HIP



INOX



COBALT CHROME



BIOLOX[®] DELTA

TECHNICAL DOCUMENTATION

12/14
HEADS RANGE



FH ORTHO heads range is available in different materials, in several diameters but with a standard 12/14 taper. All the ball heads are compatible with all the FH ORTHO stems. Regarding the liners, the ceramic liners can be match only with the ceramic heads either BIOLOX® FORTE or BIOLOX® DELTA. The TRIANON® highly crosslinked polyethylene liner can be match only with CoCr heads or with the BIOLOX® ceramic heads.


REFERENCES - NOT AVAILABLE FOR USA



STAINLESS STEEL HEADS	
231402	STAINLESS STEEL HEAD TAPER 12/14 D.22,2 M (0)
231403	STAINLESS STEEL HEAD TAPER 12/14 D.22,2 L (+2)
247852	STAINLESS STEEL HEAD TAPER 12/14 D.28 XS (-5)
231391	STAINLESS STEEL HEAD TAPER 12/14 D.28 S (-3,5)
231392	STAINLESS STEEL HEAD TAPER 12/14 D.28 M (0)
231393	STAINLESS STEEL HEAD TAPER 12/14 D.28 L (+3,5)
231394	STAINLESS STEEL HEAD TAPER 12/14 D.28 XL (+7)



COBALT CHROMIUM HEADS	
240862	COCR HEAD TAPER 12/14 D.28 S (-3,5)
240863	COCR HEAD TAPER 12/14 D.28 M (0)
240864	COCR HEAD TAPER 12/14 D.28 L (+3,5)
240865	COCR HEAD TAPER 12/14 D.28 XL (+7)
240870	COCR HEAD TAPER 12/14 D.32 S (-3,5)
240871	COCR HEAD TAPER 12/14 D.32 M (0)
240872	COCR HEAD TAPER 12/14 D.32 L (+3,5)
240873	COCR HEAD TAPER 12/14 D.32 XL (+7)



BIOLOX® DELTA HEADS	
256516	BIOLOX® DELTA HEAD TAPER 12/14 D.28 S (-3,5)
256517	BIOLOX® DELTA HEAD TAPER 12/14 D.28 M (0)
256518	BIOLOX® DELTA HEAD TAPER 12/14 D.28 L (+3,5)
256519	BIOLOX® DELTA HEAD TAPER 12/14 D.32 S (-3,5)
256520	BIOLOX® DELTA HEAD TAPER 12/14 D.32 M (0)
256521	BIOLOX® DELTA HEAD TAPER 12/14 D.32 L (+3,5)
266371	BIOLOX® DELTA HEAD TAPER 12/14 D.32 XL (+7)
256522	BIOLOX® DELTA HEAD TAPER 12/14 D.36 S (-3,5)
256523	BIOLOX® DELTA HEAD TAPER 12/14 D.36 M (0)
256524	BIOLOX® DELTA HEAD TAPER 12/14 D.36 L (+3,5)
266372	BIOLOX® DELTA HEAD TAPER 12/14 D.36 XL (+7)



DISTRIBUTEURS
DISTRIBUTORS

GRUPE
FH ORTHO™

FR, FH ORTHO SAS
3 rue de la Forêt - Zone Industrielle
BP 50009
68990 Heimsbrunn CEDEX - FRANCE
Tél. +33 (0)3 89 81 90 92
Fax : +33 (0)3 89 81 80 11
info@fhortho.com
www.fhortho.com

USA, FH ORTHOPEDICS INC.
OrthoEx
7327 E Tierra Buena Lane
Scottsdale, Arizona 85260 - USA
Phone: +1 (412) 965-0950
customerservice@fhortho-us.com
www.fhortho.com

PL, FH ORTHO POLSKA
Ul. Garbary 95/A6,
61-757 Poznan - POLSKA
Phone: +48 61 863 81 27
Fax: +48 61 863 81 28
biuro@implants24.pl
www.fhortho.com



FABRICANT
MANUFACTURER

FR, FH INDUSTRIE
6 rue Nobel, Z.I. de Kernevez
29000 QUIMPER - FRANCE
Tél. +33 (0)2 98 55 68 95
Fax : +33 (0)2 98 53 42 13
contact-fhi@fhortho.com
www.fhortho.com